

medAspis - Alert Messages

The following selection of the most common error messages appeared in the Control Panel and on the Speed Scanner.

positive responses

active This pack is able to be dispensed

successful Successfully processed

scanned The code was scanned and stored

coding alerts

unknown serial number Pack Serial Number not recognised.

Unknown Batch Number Failed to find Batch Number

Batch Number mismatch The Batch Number does not match Pack Serial Number.

Expiry Date mismatch Expiry Date does not match Batch Number.

reactivation alerts

Undo Status > 10 days
Undo timeframe was exceeded.
Undo/User mismatch
Undo attempted by different user.

Undo Status mismatch
Undo attempt doesn't match the Pack status.
Undo Status mismatch
Undo attempt doesn't match the Pack status.

double dispense alerts

inactive This pack is already inactive supplied This pack is supplied already inactive This pack is already inactive

already inactive This pack is dispensed, exported or destroyed

already inactive This pack is already inactive

sample, locked, exported, stolen, recall, etc.

destroyed This pack is already destroyed locked This pack is locked for sale exported This pack is exported sample This pack is a sample stolen This pack is marked as stolen free sample This pack is a free sample recalled This pack is recalled expired This pack is expired withdrawn This pack is withdrawn repacked This pack has been repacked

general alerts

general error A general technical exception occurred.
unknown The status of this pack is unknown
Unknown Product Code Product Code not recognised.

Double scan inside tolerance recorded.

registration alerts

invalid Arvato login
user locked
password expired
The entered ClientId or ClientSecret does not match. Invalid Arvato login
Maximum number of failed login attempts reached. User locked
The entered password is either expired or an initial password.
T&C invalid
The terms & conditions were not accepted in the allowed timeframe

double dispense recognition

dispensed here This pack has already been dispensed here

same for destroyed, locked, exported, sample, etc.

No FMD

no FMD code This code is not an FMD Code. Please check if another code exists